

Home Care: Equipment Information

Equipment Name	Must-have Inventory of Parts and Accessories	Calibration Requirement	Preventive Maintenance Schedule
Lilo			check inflation condition twice per day (early morning and late evening)
Wheelchair seat cushion			
Hospital bed and its controller			check function every week
Respirator	tubes and connectors X2		check tube and connector condition twice per day (early morning and late evening)
Oxygen Generator			check machine condition twice per day (early morning and late evening)
Oxygen Cylinder with Flow Meter Control			
Backup Battery System (or Self Power Electric Generator) for Vital System Support			
Blood Oxygen and Pulse Meter		Yes	
Phlegm(Sputum) Suction Machine	FilterX1 and Phelgm Collector X1		change filter if opaqueness occurred or too much of moisture, change collector if suction power abnormal
Stethoscope			
Inflation System Pressure Gauge Manometer			
Blood Glucose Monitoring System		Yes	
Auto Digital Blood Pressure Monitor		Yes	